

Product datasheet for **RC228601L3V**

BRD8 (NM_001164326) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | BRD8 (NM_001164326) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | BRD8 |
| Synonyms: | p120; SMAP; SMAP2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001164326 |
| ORF Size: | 2598 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC228601). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001164326.1 , NP_001157798.1 |
| RefSeq ORF: | 2601 bp |
| Locus ID: | 10902 |
| UniProt ID: | Q9H0E9 |
| Cytogenetics: | 5q31.2 |
| Protein Families: | Druggable Genome, Transcription Factors |
| MW: | 94.1 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene interacts with thyroid hormone receptor in a ligand-dependent manner and enhances thyroid hormone-dependent activation from thyroid response elements. This protein contains a bromodomain and is thought to be a nuclear receptor coactivator. Multiple alternatively spliced transcript variants that encode distinct isoforms have been identified. [provided by RefSeq, Jul 2014]