

## Product datasheet for **RC228692L3V**

### ZNF407 (NM\_001146190) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF407 (NM_001146190) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF407
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001146190
ORF Size:	4980 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228692).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001146190.1</a> , <a href="#">NP_001139662.1</a>
RefSeq ORF:	4983 bp
Locus ID:	55628
UniProt ID:	<a href="#">Q9C0G0</a>
Cytogenetics:	18q22.3
MW:	183.9 kDa
Gene Summary:	This gene encodes a zinc finger protein whose exact function is not known. It may be involved in transcriptional regulation. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]



[View online »](#)