

Product datasheet for **RC228730L4V**

SCOC (NM_001153446) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SCOC (NM_001153446) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SCOC |
| Synonyms: | HRIHFB2072; SCOCO; UNC-69 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001153446 |
| ORF Size: | 1514 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC228730). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001153446.1 , NP_001146918.1 |
| RefSeq Size: | 1876 bp |
| RefSeq ORF: | 369 bp |
| Locus ID: | 60592 |
| UniProt ID: | Q9UIL1 |
| Cytogenetics: | 4q31.1 |
| MW: | 13.9 kDa |


[View online »](#)

Gene Summary:

This gene encodes a short coiled-coiled domain-containing protein that localizes to the Golgi apparatus. The encoded protein interacts with ADP-ribosylation factor-like proteins. Pseudogenes of this gene are found on chromosomes 1 and 14. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]