

## Product datasheet for **RC228751**

### SSX2 (SSX2B) (NM\_001164417) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SSX2 (SSX2B) (NM\_001164417) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** SSX2  
**Synonyms:** CT5.2; CT5.2b; HOM-MEL-40; SSX  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC228751 representing NM\_001164417  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGAACGGAGACGACGCCTTTGCAAGGAGACCCACGGTTGGTGCTCAAATACCAGAGAAGATCCAAAAGG  
 CCTTCGATGATTTGCCAAATACTTCTCTAAGGAAGAGTGGAAAAGATGAAAGCCTCAGAGAAATCTT  
 CTATGTGTATATGAAGAGAAAGTATGAGGCTATGACTAACTAGGTTTCAAGGCCACCTCCACCTTTC  
 ATGTGTAATAAACGGGCCGAAGACTTCCAGGGGAATGATTTGGATAATGACCCTAACCGTGGGAATCAGG  
 TTGAACGTCCTCAGATGACTTTTCGGCAGGCTCCAGGGAATCTCCCGAAGATCATGCCAAGAAGCCAGC  
 AGAGGAAGGAAATGATTCGGAGGAAGTCCAGAAAGCATCTGGCCACAAAATGATGGGAAAGAGCTGTGC  
 CCCCCGGGAAAACCAACTACCTCTGAGAAGATTACGAGAGATCTGGACCCAAAAGGGGGGAACATGCCT  
 GGACCCACAGACTGCGTGAGAGAAAACAGCTGGTGATTTATGAAGAGATCAGCGACCCTGAGGAAGATGA  
 CGAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC228751 representing NM\_001164417  
 Red=Cloning site Green=Tags(s)

MNGDDAFARRPTVGAQIPEKIQAFFDDIAKYFSKEEWEKMKASEKIFVYVMKRKYEAMTKLGFKATLPPF  
 MCNKRAEDFQGNLDNDPNRGNQVERPQMTFGRLQGISPKIMPKKPAEEGNDSEEVPEASGPQNDGKELC  
 PPGKPTTSEKIHRSRGPKRGEHAWTHRLRERKQLVIYEEISDPEEDDE

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV




[View online »](#)

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja1446\\_b09.zip](https://cdn.origene.com/chromatograms/ja1446_b09.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001164417

**ORF Size:** 564 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001164417.3](#)

**RefSeq ORF:** 567 bp

**Locus ID:** 727837


**UniProt ID:** [Q16385](#)

**Cytogenetics:** Xp11.22


**MW:** 21.4 kDa

**Gene Summary:** The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. This gene, and also the SSX1 and SSX4 family members, have been involved in t(X;18)(p11.2;q11.2) translocations that are characteristically found in all synovial sarcomas. This translocation results in the fusion of the synovial sarcoma translocation gene on chromosome 18 to one of the SSX genes on chromosome X. The encoded hybrid proteins are likely responsible for transforming activity. Alternative splicing of this gene results in multiple transcript variants. This gene also has an identical duplicate, GeneID: 6757, located about 45 kb upstream in the opposite orientation on chromosome X. [provided by RefSeq, Jul 2013]


**Product images:**



Circular map for RC228751



Western blot validation of overexpression lysate (Cat# [LY431779]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC228751 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified SSX2B protein (Cat# [TP328751]). The protein was produced from HEK293T cells transfected with SSX2B cDNA clone (Cat# RC228751) using MegaTran 2.0 (Cat# [TT210002]).