

Product datasheet for **RC228763A1V**

Human ABHD14B (NM_001146314) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ABHD14B (NM_001146314) AAV Particle
Tag: Myc-DDK
Symbol: ABHD14B
Synonyms: CIB; HEL-S-299
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC228763 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGCAAGCGTGGAGCAGCGCAGGGCACCATCCAGGTGCAGGGCCAGGCCCTCTTCTTCCGAGAGG
CCCTGCCCCGCGAGTGGGAGGCTCGCTTCTCTGTACTGCTGCTGCATGGTATTCGCTTCTCTCCGAGAC
CTGGCAGAACCTGGGTACACTGCACAGGCTGGCCAGGCTGGCTACCGGGCTGTGGCCATTGACCTGCCA
GGTCTGGGGCACTCCAAGGAAGCAGCAGCCCCGCCCCATTGGGGAGCTGGCCCTGGCAGCTTCCTGG
CGGCTGTGGTGGATGCCCTGGAGCTGGGCCCCCGGTTGTGATCAGTCCATCACTGAGTGGCATGTACTC
CCTGCCCTTCTCACGGCCCCGCTGCCAGCTCCCGGGCTTTGTGCCAGTGGCCCCCATCTGCAGTGC
AAAATCAATGCTGCCAATATGCCAGTGTGAAGACTCCAGCTCTGATTGTATATGGAGACCAGGACCCCA
TGGGTGAGACCAGCTTTGAGCACCTGAAGCAGCTGCCCAACCACCGGGTGTGATCATGAAGGGGCGGG
GCACCCCTGTTACCTGGACAAACCAGAGGAGTGGCATACAGGGCTGCTGGACTTCCTGCAGGGGCTCCAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228763 protein sequence
Red=Cloning site Green=Tags(s)

MAASVEQREGTIQVQGQALFFREALPGSGQARFSVLLHGIRESSETWQNLGTLHRLAQAGYRAVAIDL
GLGHSKEAAAPAPIGELAPGSFLAAVDALELGPPVVISPSLSGMYSLPFLTAPGSQPLPGFVPVAPICTD
KINAANYASVKTPALIVYGDQDPMGQTSFEHLKQLPNHRVLIMKGAGHPCYLDKPEEWHGLLDLFLQGLQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001146314
ORF Size:	630 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001146314.1</u> , <u>NP_001139786.1</u>
RefSeq Size:	2253 bp
RefSeq ORF:	633 bp
Locus ID:	84836
UniProt ID:	<u>Q96IU4</u>
Cytogenetics:	3p21.2
MW:	22.3 kDa