

Product datasheet for **RC228763L4V**

ABHD14B (NM_001146314) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ABHD14B (NM_001146314) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABHD14B
Synonyms:	CIB; HEL-S-299
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001146314
ORF Size:	630 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228763).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001146314.1 , NP_001139786.1
RefSeq Size:	2253 bp
RefSeq ORF:	633 bp
Locus ID:	84836
UniProt ID:	Q96IU4
Cytogenetics:	3p21.2
MW:	22.3 kDa



[View online »](#)

Gene Summary: Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.
[UniProtKB/Swiss-Prot Function]