

## Product datasheet for **RC228770L3V**

### Secretogranin 3 (SCG3) (NM\_001165257) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Secretogranin 3 (SCG3) (NM_001165257) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCG3
Synonyms:	SGIII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001165257
ORF Size:	1407 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228770).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001165257.1</a> , <a href="#">NP_001158729.1</a>
RefSeq Size:	3313 bp
RefSeq ORF:	711 bp
Locus ID:	29106
UniProt ID:	<a href="#">Q8WXD2</a>
Cytogenetics:	15q21.2
Protein Families:	Druggable Genome, Secreted Protein
MW:	53 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Granins may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]