

## Product datasheet for **RC228808L4V**

### **NAT1 (NM\_001160170) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NAT1 (NM_001160170) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NAT1
Synonyms:	AAC1; MNAT; NAT-1; NATI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001160170
ORF Size:	870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228808).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001160170.1</a>
RefSeq Size:	2355 bp
RefSeq ORF:	873 bp
Locus ID:	9
UniProt ID:	<a href="#">P18440</a>
Cytogenetics:	8p22
Protein Pathways:	Caffeine metabolism, Drug metabolism - other enzymes, Metabolic pathways
MW:	33.9 kDa



[View online »](#)

**Gene Summary:**

This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]