

## Product datasheet for **RC228809**

### **NAT1 (NM\_001160171) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NAT1 (NM_001160171) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NAT1
Synonyms:	AAC1; MNAT; NAT-1; NATI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC228809 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGACATTGAAGCATATCTTGAAAGAATTGGCTATAAGAAGTCTAGGAACAAATTGGACTTGGAAACAT  
TAACTGACATTCTTCAACACCAGATCCGAGCTGTTCCCTTTGAGAACCTTAACATCCATTGTGGGGATGC  
CATGGACTTAGGCTTAGAGGCCATTTTTGATCAAGTTGTGAGAAGAAATCGGGTGGATGGTGTCTCCAG  
GTCAATCATCTTCTGTACTGGGCTCTGACCACTATTGGTTTTGAGACCACGATGTTGGGAGGGTATGTTT  
ACAGCACTCCAGCCAAAAATACAGCACTGGCATGATTCACCTTCTCCTGCAGGTGACCATTGATGGCAG  
GAACTACATTGTCGATGCTGGGTTTGGACGCTCATACCAGATGTGGCAGCCTCTGGAGTTAATTTCTGGG  
AAGGATCAGCCTCAGGTGCCTTGTGCTTCCGTTTGACGGAAGAGAATGGATTCTGGTATCTAGACCAAA  
TCAGAAGGGAACAGTACATTCCAATGAAGAATTTCTTCATTCTGATCTCCTAGAAGACAGCAAATACCG  
AAAAATCTACTCCTTACTCTTAAGCCTCGAACAATTGAAGATTTTGGTCTATGAATACATACCTGCAG  
ACATCTCCATCATCTGTGTTTACTAGTAAATCATTTTGTCTTGCAGACCCAGATGGGGTTCACTGTT  
TGGTGGGCTTCAACCTACCCATAGGAGATTCAATTATAAGGACAATACAGATCTAATAGAGTTCAAGAC  
TCTGAGTGAGGAAGAAATAGAAAAAGTCTGAAAAATATATTTAATATTTCCCTTGCAGAGAAAGCTTGTG  
CCCAAACATGGTGATAGATTTTTACTATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC228809 protein sequence  
Red=Cloning site Green=Tags(s)

MDIEAYLERIGYKKS RNKLDLETLDILQHQIRAVPFENLNIHCGDAMD LGLEAIFDQVRRNRGGWCLQ  
 VNHLLYWALTTIGFETTMLGGYVYSTPAKKYSTGMIHLLQVTIDGRNYIVDAGFGRSYQMWQPLELISG  
 KDQPQVPCVFRLTEENGFWYLDQIRREQYIPNEEFLHSDLLLEDSKYRKIYSFTLKPRTIEDFESMNTYLQ  
 TSPSSVFTSKSFCSLQTPDGVHCLVGF TLTHRRFNYKDNTDLIEFKTLSEEEIEKVLKNI FNISLQRKLV  
 PKHGDRFFTI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8075\\_f10.zip](https://cdn.origene.com/chromatograms/mk8075_f10.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001160171

**ORF Size:** 870 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001160171.4](#)

**RefSeq Size:** 2205 bp

**RefSeq ORF:** 873 bp

**Locus ID:** 9

**UniProt ID:** [P18440](#)


**Cytogenetics:** 8p22

**Protein Pathways:** Caffeine metabolism, Drug metabolism - other enzymes, Metabolic pathways


**MW:** 33.9 kDa

**Gene Summary:** This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]

**Product images:**



Circular map for RC228809



Western blot validation of overexpression lysate (Cat# [LY431840]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC228812] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).