

## Product datasheet for **RC228818L3V**

### NEK6 (NM\_001166168) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NEK6 (NM_001166168) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEK6
Synonyms:	SID6-1512
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001166168
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228818).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001166168.1</a> , <a href="#">NP_001159640.1</a>
RefSeq Size:	2587 bp
RefSeq ORF:	942 bp
Locus ID:	10783
UniProt ID:	<a href="#">Q9HC98</a>
Cytogenetics:	9q33.3
Protein Families:	Druggable Genome, Protein Kinase
MW:	35.7 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a kinase required for progression through the metaphase portion of mitosis. Inhibition of the encoded protein can lead to apoptosis. This protein also can enhance tumorigenesis by suppressing tumor cell senescence. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]