

Product datasheet for **RC228828A1V**

Human GPR17 (NM_001161417) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human GPR17 (NM_001161417) AAV Particle
Tag:	Myc-DDK
Symbol:	GPR17
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC228828 representing NM_001161417 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAATGGCCTTGAAGTGGCTCCCCAGGTCTGATCACCAACTTCTCCCTGGCCACGGCAGAGCAATGTG
GCCAGGAGACGCCACTGGAGAACATGCTGTTGCCTCCTTCTACCTTCTAGATTTATCCTGGCTTTAGT
TGGCAATACCCTGGCTCTGTGGTTTTTCATCCGAGACCACAAGTCCGGGACCCCGGCCAACGTGTTCTGT
ATGCATCTGGCCGTGGCCGACTTGTCTGCGTGTGGTCTGCCACCCGCTGGTCTACCACTTCTCTG
GGAACCACTGGCCATTTGGGAAATCGCATGCCGTCTCACCAGGCTTCTCTTCTACCTCAACATGTACGC
CAGCATCTACTTCCTCACCTGCATCAGCGCCGACCGTTTCTGGCCATTGTGCACCCGGTCAAGTCCCTC
AAGCTCCGCAGGCCCTCTACGCACACCTGGCCTGTGCCTTCTGTGGGTGGTGGTGGCTGTGGCCATGG
CCCCGCTGCTGGTGAGCCACAGACCGTGCAGACCAACCACACGGTGGTCTGCCTGCAGCTGTACCGGGA
GAAGGCCTCCACCATGCCCTGGTGTCCCTGGCAGTGGCCTTCACCTTCCCGTTCATCACCACGGTCACC
TGCTACCTGCTGATCATCCGACGCTGCGGCAGGGCCTGCGTGTGGAGAAGCGCCTCAAGACCAAGGCAG
TGCGCATGATCGCCATAGTGTGGCCATCTTCTGGTCTGCTTCGTGCCCTACCACGTCAACCGCTCCGT
CTACGTGCTGCACTACCGCAGCCATGGGGCCTCCTGCGCCACCCAGCGCATCCTGGCCCTGGCAAACCGC
ATCACCTCCTGCCTCACCAGCCTCAACGGGGCACTCGACCCCATCATGTATTTCTCGTGGCTGAGAAGT
TCCGCCACGCCCTGTGCAACTTGCTCTGTGGCAAAGGCTCAAGGGCCCGCCCCCAGCTTCGAAGGGAA
AACCAACGAGAGCTCGCTGAGTGCCAAGTCAGAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC228828 representing NM_001161417 Red=Cloning site Green=Tags(s) MNGLEVAPPGLITNFSLATAEQCGQETPLENMLFASFYLLDFILALVGNTLALWFFIRDHKSGTPANVFL MHLAVADLSCVLVLPTRLVYHFSGNHWPFGIACRLTGFLFYLNMYASIYFLTCTISADRFLAIVHPVKSL KLRRPLYAHLACAFLWVVAVAMAPLLVSPQTVQTNHTVVCLQLYREKASHHALVSLAVAFTFPFITTVT CYLLIIRSLRQGLRVEKRLKTKAVRMIAIVLAIFLVCFVPYHVNRSVYVLHYRSHGASCATQRILALANR ITSCLTSLNGALDPIMYFFVAEKFRHALCNLLCGKRLKGPPPSFEGKTNESSLSAKSEL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001161417
ORF Size:	1017 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001161417.1</u>
RefSeq Size:	2315 bp
RefSeq ORF:	1020 bp
Locus ID:	2840
UniProt ID:	<u>Q13304</u>
Cytogenetics:	2q14.3
MW:	37.8 kDa