

Product datasheet for **RC228834L2V**

P2RY5 (LPAR6) (NM_001162497) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | P2RY5 (LPAR6) (NM_001162497) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | LPAR6 |
| Synonyms: | ARWH1; HYPT8; LAH3; LPA-6; P2RY5; P2Y5 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001162497 |
| ORF Size: | 1032 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC228834). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001162497.1 |
| RefSeq Size: | 2561 bp |
| RefSeq ORF: | 1035 bp |
| Locus ID: | 10161 |
| UniProt ID: | P43657 |
| Cytogenetics: | 13q14.2 |
| Protein Families: | Druggable Genome, GPCR, Transmembrane |
| Protein Pathways: | Neuroactive ligand-receptor interaction |



[View online »](#)

MW: 39.4 kDa

Gene Summary: The protein encoded by this gene belongs to the family of G-protein coupled receptors, that are preferentially activated by adenosine and uridine nucleotides. This gene aligns with an internal intron of the retinoblastoma susceptibility gene in the reverse orientation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2009]