

Product datasheet for **RC228850L4V**

GPR85 (NM_001146266) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | GPR85 (NM_001146266) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | GPR85 |
| Synonyms: | SREB; SREB2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001146266 |
| ORF Size: | 1110 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC228850). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001146266.1 , NP_001139738.1 |
| RefSeq Size: | 4841 bp |
| RefSeq ORF: | 1113 bp |
| Locus ID: | 54329 |
| UniProt ID: | P60893 |
| Cytogenetics: | 7q31.1 |
| Protein Families: | Druggable Genome, GPCR, Transmembrane |
| MW: | 42 kDa |



[View online »](#)

Gene Summary:

Members of the G protein-coupled receptor (GPCR) family, such as GPR85, have a similar structure characterized by 7 transmembrane domains. Activation of GPCRs by extracellular stimuli, such as neurotransmitters, hormones, or light, induces an intracellular signaling cascade mediated by heterotrimeric GTP-binding proteins, or G proteins (Matsumoto et al., 2000 [PubMed 10833454]).[supplied by OMIM, Aug 2008]