

Product datasheet for **RC22853A1V**

Human LRTM2 (NM_001163925) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LRTM2 (NM_001163925) AAV Particle
Tag:	Myc-DDK
Symbol:	LRTM2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC22853 representing NM_001163925 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTGGCGCCGGGCAGCAGCCCTGGGCAGAGGGGCAGGCTCGCCCTGCAGTGGAGGCAAGTCTCCTGGA
TCACCTGCTGGATCGCCCTGTATGCTGTGGAGGCCCTCCACCTGCCCTTCTCCTGCAAGTGTGACAG
CCGCAGCCTGGAGGTGGACTGCAGTGGCCTTGGCCTCACCACGGTGCCCCAGACGTGCCCGCAGCCACC
CGAACCCCTCTTGCTCTTGAACAATAAGCTGAGTGCCCTGCCAAGCTGGGCTTTCGCCAACCTCTCCAGCC
TGCAGCGGTTGGACCTGTCCAACAACCTCCTGGACCGGCTGCCCGCTCCATTTTCGGGGACCTGACGAA
TCTGACTGAGCTTCAGCTGCGCAATAACAGCATCAGGACCCTGGACAGGACCTGCTGCGGCACTCGCCG
CTGCTCCGCCACCTGGACCTGTCCATCAACGGCCTGGCCAGTTGCCCCCTGGTCTTTTCGACGGGCTCC
TGGCTCTGCGCTCCCTCTCGCTTCCGCTCCAACCGTCTGCAGAATCTGGACCGGCTGACATTTGAACCCCT
AGCAAACCTGCAGCTGCTGCAGGTGGGGATAACCCCTGGGAGTGTGACTGTAACCTGCGTGAAGTTCAA
CACTGGATGGAGTGGTTCTCCTACCGAGGGGGACGCTTGGACCAGCTTGCTGCACCTGCCCAAGGAGC
TGAGGGGAAGGACATGCGGATGGTCCCATGGAGATGTTCAACTACTGCTCCAGCTGGAGGACGAGAA
TAGCTCAGCTGGGCTGGATATTCCTGGGCCACCCTGCACCAAGGCCAGTCCAGAGCCTGCTAAGCCCAAG
CCCGGGGCTGAGCCGGAGCCGGAGCCAGCACAGCCTGCCACAGAAGCAGAGGACCCGCCGGCGAGCG
TGAGGCGAGCCATGGGCACGGTGATCATTGCAGGGGTCTGTGCGGCGTCTGTCATCATGATGGTGGT
GGCCGCTGCCTATGGCTGCATCTACGCCTCCCTCATGGCCAAGTACCACCGGGAGCTCAAAAAGCGCCAG
CCCCTGATGGGGACCCCGAGGGCGAGCACGAGGACCAGAAGCAGATCTTCTGTGGCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC228853 representing NM_001163925 Red=Cloning site Green=Tags(s) MLAPGSSPGQRGRLALQWRQVSWITCWIALYAVEALPTCPFSCCKDSRSLEVDCSGLGLTTVPPDVPAAT RTLLLLNNKLSALPSWAFANLSSLQRLDLNNFLDRLPRSIFGDLTNLTELQLRNNSIRTLDRLLRHSP LLRHLDLSINGLAQLPPGLFDGLLALRSLSLRSNRLQNLDRLTFEPLANLQLQVGNPWECDCNREFK HMEWFSYRGGRLDQLACTLPKELRGKDMRMVPMEMFNYSQLEDENSSAGLDIPGPPCTKASPEPAKPK PGAEPPEPSTACPQKQRHRPASVRRAMGTVIIAGVVCVVCIIMMVAAAYGCIYASLMAKYHRELKKRQ PLMGDPEGEHEDQKQISSVA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001163925
ORF Size:	1110 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001163925.1, NP_001157397.1</u>
RefSeq ORF:	1113 bp
Locus ID:	654429
UniProt ID:	<u>Q8N967</u>
Cytogenetics:	12p13.33
MW:	41 kDa