

Product datasheet for **RC22859A1V**

Human PRKAR1B (NM_001164758) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PRKAR1B (NM_001164758) AAV Particle
Tag:	Myc-DDK
Symbol:	PRKAR1B
Synonyms:	PRKAR1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC22859 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCTCCCCGCCGCTGCCCTCGGAGGAGGACGAGAGCCTGAAGGGCTGTGAGCTGTACGTGCAGC
TGCACGGGATCCAGCAGGTCTCAAAGACTGTATCGTCCACCTCTGCATCTCCAAGCCCCAACGCCCAT
GAAGTTCTCCGGGAGCACTTCGAGAAGCTGGAGAAGGAAGAAACAGGCAGATTTTGGCGCGGCAAAAG
TCAAATCACAGTCGGACTCCCATGATGAGGAGGTGTGCGCCACCCCCCGAACCTGTGGTGAAGGCC
GCCGCCGGCGAGGAGGCGTGAGTGCCGAGGTGTACACCGAGGAGGACGCCGTGTCCTACGTCAGGAAGGT
GATTCCTCAAGGACTACAAAACCATGACTGCGCTGGCCAAAGGCCATCTCCAAGAAGCTGCTCTTCGCTCAC
CTGGATGACAACGAGAGGAGTGACATATTCGATGCCATGTTCCCTGTCACTCACATCGCTGGGAGACTG
TTATACAGCAAGGGAATGAAGGAGACAATTCTATGTCGTTGATCAAGGGGAAGTGGATGTGTACGTGAA
CGGAGAGTGGGTGACCAACATCAGCGAGGGAGGAGCTTCGGGGAGCTGGCGCTCATCTACGGCACCCCC
AGGGCTGCGACCGTGAAAGCCAAGACGGACCTCAAGCTCTGGGGGATCGACCGGGACAGCTACCGCGCA
TCCTTATGGGCAGCACGCTGAGGAAACGAAGATGTACGAGGAGTTCTCAGCAAGGTCTCCATCTAGA
GTCCCTGGAGAAGTGGGAGCGTCTGACCGTGGCGGATGCGCTGGAGCCGTCAGTTTGAAGATGGAGAG
AAAATCGTGGTCCAGGGAGAGCCTGGGACGACTTTTACATCATCACGAGGGCACCGCGTCCGTGCTGC
AGCGCCCGTCCCCAATGAGGAGTACGTGGAGGTGGGGCGCTGGGACCTCTGACTACTTCGGGAGAT
TGCACTGCTGCTGAACCGGCCCGGGCGGCACTGTGCTGGCCCGGGGGCCCTCAAGTGTGTGAAGCTG
GACCGGCCCGCTTCGAGCGTGTGCTGGGGCCCTGCTCTGAGATCCTCAAGAGGAACATTCAGCGTTACA
ACAGCTTCATCTCCCTCACCGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC228859 protein sequence Red=Cloning site Green=Tags(s) MASPPACPSEEDESLKGCELYVQLHGIQQVLKDCIVHLCISKPERPMKFLREHFEEKLEKEENRQILARQK SNSQSDSHDEEVSPTPNPVVKARRRRGGVSAEVYTEEDAVSYVRKVIPKDYKTMTALAKAISKNVLF AH LDDNERSDIFDAMFPVTHIAGETVIQQGNEGDNFYVVDQGEVDVYVNGEWVTNISEGGSFGELALIYGTP RAATVKAKTDLKLWGIDRDSYRRILMGSTLRKRKMYEEFLSKVSILESEKWERLTVADALEPVQFEDGE KIVVQGEPGDDFYIITEGTASVLQRRSPNEEYVEVGRLGPSDYFGEIALLLNRPRAATTVARGPLKCVKL DRPRFERVLGPCSEILKRNIQRYNSFISLTV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001164758
ORF Size:	1143 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001164758.1</u> , <u>NP_001158230.1</u>
RefSeq Size:	2488 bp
RefSeq ORF:	1146 bp
Locus ID:	5575
UniProt ID:	<u>P31321</u>
Cytogenetics:	7p22.3
MW:	43.1 kDa