

Product datasheet for **RC228876L3V**

LENG4 (MBOAT7) (NM_001146056) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LENG4 (MBOAT7) (NM_001146056) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LENG4
Synonyms:	BB1; hMBOA-7; LENG4; LPIAT; LPLAT; LRC4; MBOA7; MRT57; OACT7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001146056
ORF Size:	1197 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228876).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001146056.1
RefSeq ORF:	1200 bp
Locus ID:	79143
UniProt ID:	Q96N66
Cytogenetics:	19q13.42
Protein Families:	Transmembrane
MW:	44.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the membrane-bound O-acyltransferases family of integral membrane proteins that have acyltransferase activity. The encoded protein is a lysophosphatidylinositol acyltransferase that has specificity for arachidonoyl-CoA as an acyl donor. This protein is involved in the reacylation of phospholipids as part of the phospholipid remodeling pathway known as the Land cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2009]