

## Product datasheet for **RC228906L4V**

### TCF7L2 (NM\_001146284) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | TCF7L2 (NM_001146284) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | TCF7L2   |
| Synonyms:                 | TCF-4; TCF4  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001146284   |
| ORF Size:                 | 1365 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC228906).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001146284.1</a>   |
| RefSeq ORF:               | 1368 bp  |
| Locus ID:                 | 6934   |
| UniProt ID:               | <a href="#">Q9NQB0</a>   |
| Cytogenetics:             | 10q25.2-q25.3  |
| Protein Families:         | Druggable Genome, Transcription Factors  |
| Protein Pathways:         | Acute myeloid leukemia, Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Basal cell carcinoma, Colorectal cancer, Endometrial cancer, Melanogenesis, Pathways in cancer, Prostate cancer, Thyroid cancer, Wnt signaling pathway  |



[View online »](#)

**MW:** 50.4 kDa

**Gene Summary:** This gene encodes a high mobility group (HMG) box-containing transcription factor that plays a key role in the Wnt signaling pathway. The protein has been implicated in blood glucose homeostasis. Genetic variants of this gene are associated with increased risk of type 2 diabetes. Several transcript variants encoding multiple different isoforms have been found for this gene.[provided by RefSeq, Oct 2010]