

Product datasheet for RC228921L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Calcitonin receptor (CALCR) (NM_001164738) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Calcitonin receptor (CALCR) (NM_001164738) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Calcitonin receptor

Synonyms: CRT; CT-R; CTR1

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC228921).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_001164738

ORF Size: 1422 bp



Calcitonin receptor (CALCR) (NM_001164738) Human Tagged Lenti ORF Clone - RC228921L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001164738.1</u>, <u>NP 001158210.1</u>

7q21.3

RefSeq ORF: 1425 bp

Locus ID: 799

Cytogenetics:

UniProt ID: P30988

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Neuroactive ligand-receptor interaction

MW: 55.34 kDa


Gene Summary: This gene encodes a high affinity receptor for the peptide hormone calcitonin and belongs to

a subfamily of seven transmembrane-spanning G protein-coupled receptors. The encoded protein is involved in maintaining calcium homeostasis and in regulating osteoclast-mediated bone resorption. Polymorphisms in this gene have been associated with variations in bone mineral density and onset of osteoporosis. Alternate splicing results in multiple transcript

variants. [provided by RefSeq, Sep 2009]



Product images:



Circular map for RC228921L4