

Product datasheet for RC228924L1V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

AMHR2 (NM_001164690) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: AMHR2 (NM_001164690) Human Tagged ORF Clone Lentiviral Particle

Symbol: AMHR2

Synonyms: AMHR; MISR2; MISRII; MRII

Mammalian Cell

Selection:

None

Vector:

pLenti-C-Myc-DDK (PS100064)

Tag: Myc-DDK

ACCN: NM_001164690

ORF Size: 1434 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC228924).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 001164690.1

RefSeq Size:1859 bpRefSeq ORF:1437 bpLocus ID:269

 UniProt ID:
 Q16671

 Cytogenetics:
 12q13.13

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction, TGF-beta signaling pathway





ORIGENE

MW: 52.2 kDa

Gene Summary:

This gene encodes the receptor for the anti-Mullerian hormone (AMH) which, in addition to testosterone, results in male sex differentiation. AMH and testosterone are produced in the testes by different cells and have different effects. Testosterone promotes the development of male genitalia while the binding of AMH to the encoded receptor prevents the development of the mullerian ducts into uterus and Fallopian tubes. Mutations in this gene are associated with persistent Mullerian duct syndrome type II. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Sep 2009]