

Product datasheet for **RC228971L3V**

BEGAIN (NM_001159531) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BEGAIN (NM_001159531) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BEGAIN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001159531
ORF Size:	1779 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228971).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001159531.1 , NP_001153003.1
RefSeq Size:	2763 bp
RefSeq ORF:	1782 bp
Locus ID:	57596
UniProt ID:	Q9BUH8
Cytogenetics:	14q32.2
Protein Families:	Druggable Genome
MW:	64.8 kDa
Gene Summary:	May sustain the structure of the postsynaptic density (PSD).[UniProtKB/Swiss-Prot Function]



[View online »](#)