

## Product datasheet for **RC228980L4V**

### RPAP3 (NM\_001146075) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	RPAP3 (NM_001146075) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RPAP3
Synonyms:	hSpagh; Tah1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001146075
ORF Size:	1893 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228980).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001146075.1</a> , <a href="#">NP_001139547.1</a>
RefSeq ORF:	1896 bp
Locus ID:	79657
UniProt ID:	<a href="#">Q9H6T3</a>
Cytogenetics:	12q13.11
Protein Families:	Protease
MW:	71.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an RNA polymerase II-associated protein. The encoded protein may function in transcriptional regulation and may also regulate apoptosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]