

## Product datasheet for **RC229034L4V**

### ALIX (PDCD6IP) (NM\_001162429) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ALIX (PDCD6IP) (NM_001162429) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PDCD6IP
Synonyms:	AIP1; ALIX; DRIP4; HP95
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001162429
ORF Size:	2619 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229034).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001162429.2</a> , <a href="#">NP_001155901.1</a>
RefSeq Size:	5987 bp
RefSeq ORF:	2622 bp
Locus ID:	10015
UniProt ID:	<a href="#">Q8WUM4</a>
Cytogenetics:	3p22.3
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis



[View online »](#)

MW: 96.8 kDa

**Gene Summary:** This gene encodes a protein that functions within the ESCRT pathway in the abscission stage of cytokinesis, in intraluminal endosomal vesicle formation, and in enveloped virus budding. Studies using mouse cells have shown that overexpression of this protein can block apoptosis. In addition, the product of this gene binds to the product of the PDCD6 gene, a protein required for apoptosis, in a calcium-dependent manner. This gene product also binds to endophilins, proteins that regulate membrane shape during endocytosis. Overexpression of this gene product and endophilins results in cytoplasmic vacuolization, which may be partly responsible for the protection against cell death. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. Related pseudogenes have been identified on chromosome 15. [provided by RefSeq, Jan 2012]