

Product datasheet for **RC229035L3V**

SAPS3 (PPP6R3) (NM_001164161) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SAPS3 (PPP6R3) (NM_001164161) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SAPS3 |
| Synonyms: | C11orf23; PP6R3; SAP190; SAPL; SAPLa; SAPS3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001164161 |
| ORF Size: | 2619 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC229035). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001164161.1 |
| RefSeq ORF: | 2622 bp |
| Locus ID: | 55291 |
| UniProt ID: | Q5H9R7 |
| Cytogenetics: | 11q13.2 |
| MW: | 97.5 kDa |



[View online »](#)

Gene Summary:

Protein phosphatase regulatory subunits, such as SAPS3, modulate the activity of protein phosphatase catalytic subunits by restricting substrate specificity, recruiting substrates, and determining the intracellular localization of the holoenzyme. SAPS3 is a regulatory subunit for the protein phosphatase-6 catalytic subunit (PPP6C; MIM 612725) (Stefansson and Brautigan, 2006 [PubMed 16769727]).[supplied by OMIM, Nov 2010]