

Product datasheet for **RC229073L4V**

NTE (PNPLA6) (NM_001166113) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NTE (PNPLA6) (NM_001166113) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NTE
Synonyms:	BNHS; iPLA2delta; LNMS; NTE; NTEMND; OMCS; SPG39; sws
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001166113
ORF Size:	3981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229073).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001166113.1 , NP_001159585.1
RefSeq ORF:	3984 bp
Locus ID:	10908
UniProt ID:	Q8IY17
Cytogenetics:	19p13.2
Protein Families:	Transmembrane
MW:	146 kDa



[View online »](#)

Gene Summary:

This gene encodes a phospholipase that deacetylates intracellular phosphatidylcholine to produce glycerophosphocholine. It is thought to function in neurite outgrowth and process elongation during neuronal differentiation. The protein is anchored to the cytoplasmic face of the endoplasmic reticulum in both neurons and non-neuronal cells. Mutations in this gene result in autosomal recessive spastic paraplegia, and the protein is the target for neurodegeneration induced by organophosphorus compounds and chemical warfare agents. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]