

Product datasheet for **RC229118A1V**

Human ZNF268 (NM_152943) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ZNF268 (NM_152943) AAV Particle
Tag: Myc-DDK
Symbol: ZNF268
Synonyms: HZF3
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC229118 representing NM_152943
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCACCAGGGTCCGACAGCTTCTATTTGGGTCCACCTCTCCAAGAACGAAACAGTTCATGGGATA
GGATCAGAAAGCTCCAAGGTCAGGAATCCATCTTGGGCCAAGGGACTCCTGGTCTGCAACCTCTCCCTGG
AACACCCAGGCAGAAGCAGAAGAGTCGAGAAATAGAGAAAGTCTAGAGTGGCTGTTTATTTCCCAAGAG
CAGCCAAAAATCACCAAGTCTTGGGGACCTTTGTCATTATGATGTGTTTGTGGATTTTACCTGGGAGG
AGTGGCAGCTGCTAGACCCAGCACAGAAGTGCCTGTACAGGAGTGTGATGTTGGAGAAGTATAGCAACCT
GGTGTCCCTAGGGTACCAACACACCAAACTGATATCATCTTCAAGTTGGAACAAGGAGAAGAGCTGTGT
ATGGTGCAGGCCCAAGTTCCAAATCAGACCTGTCCAATTTGAAGGCTGAAAAGTCCAAAGCCAAGGTGC
TGGCAGGTTTGGTGTCTGGTGAAGGCTGCTCTGTGCTTCCAAGATGACGCCTTGTGCTGCATCCTCTG
GAGACACAGTCTGAAAAAT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229118 representing NM_152943
Red=Cloning site **Green**=Tags(s)

MATRVRTASIVPPLQERNSSWDRIRKLQGQESILGQGTPLQPLPGTPRQKQKSRRIEKVLEWLFISQE
QPKITKSWGPLSFMDVFDFTWEEWQLLDPAQKCLYRSVMLENYSNLVSLGYQHTKPDIIIFKLEQGEELC
MVQAQVPNQTCPIKAGKSKAKVLAVLSGEGLLCASKMTPCCILWRHLEN

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_152943
ORF Size:	579 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152943.2</u>
RefSeq ORF:	582 bp
Locus ID:	10795
UniProt ID:	<u>Q14587</u>
Cytogenetics:	12q24.33
MW:	21.8 kDa