

Product datasheet for **RC229138L4V**

FCAMR (NM_032029) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FCAMR (NM_032029) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FCAMR
Synonyms:	CD351; FCA/MR; FKSG87
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032029
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229138).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032029.4 , NP_114418.2
RefSeq Size:	1981 bp
RefSeq ORF:	798 bp
Locus ID:	83953
UniProt ID:	Q8WWV6
Cytogenetics:	1q32.1
Protein Families:	Druggable Genome, Transmembrane
MW:	29.9 kDa



[View online »](#)

Gene Summary:

Functions as a receptor for the Fc fragment of IgA and IgM. Binds IgA and IgM with high affinity and mediates their endocytosis. May function in the immune response to microbes mediated by IgA and IgM.[UniProtKB/Swiss-Prot Function]