

Product datasheet for **RC229147L3V**

PPM1M (NM_144641) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPM1M (NM_144641) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPM1M
Synonyms:	PP2C-eta; PP2CE; PP2Ceta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144641
ORF Size:	894 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229147).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144641.2 , NP_653242.2
RefSeq ORF:	1380 bp
Locus ID:	132160
UniProt ID:	Q96MI6
Cytogenetics:	3p21.2
Domains:	PP2C
Protein Families:	Druggable Genome, Phosphatase
MW:	33.3 kDa



[View online »](#)