

## Product datasheet for **RC229215L4V**

### PNCK (NM\_001039582) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PNCK (NM_001039582) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PNCK
Synonyms:	BSTK3; CaMK1b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001039582
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229215).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039582.3</a> , <a href="#">NP_001034671.3</a>
RefSeq ORF:	1281 bp
Locus ID:	139728
UniProt ID:	<a href="#">Q6P2M8</a>
Cytogenetics:	Xq28
Protein Families:	Druggable Genome, Protein Kinase
MW:	46.8 kDa



[View online »](#)

**Gene Summary:**

PNCK is a member of the calcium/calmodulin-dependent protein kinase family of protein serine/threonine kinases (see CAMK1; MIM 604998) (Gardner et al., 2000 [PubMed 10673339]).[supplied by OMIM, Mar 2008]