

Product datasheet for **RC229221A1V**

Human Caspase 5 (CASP5) (NM_004347) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Caspase 5 (CASP5) (NM_004347) AAV Particle
Tag:	Myc-DDK
Symbol:	Caspase 5
Synonyms:	ICE(rel)III; ICEREL-III; ICH-3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC229221 representing NM_004347 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTGAAGACAGTGGCAAAAAAAAAAGGCGTAAGAATTTGAAGCTATGTTCAAAGGTATCCTTCAGAGTGGATTGGATAA
CTTCGTGATAAACACCATGCTAAAGAACACGTGGCTGGACAAACATCTATCCAGAC
CCTAGTACCTAATACGGATCAAAAGTCGACCAAGTGTAAAAAGACAACCACAAAAAACAGTTAAG
ATGTTGGAATACCTGGCAAGATGTTCTTCATGGTGTTTTAATTATTTGGCAAAACACGATGTTCTGA
CATTGAAGGAAGAGGAAAAAGAAAAATATTATGATACCAAAATTGAAGACAAGGCCCTGATCTTGGTAGA
CTCTTTTGCAGAAAGAAATCGCGTGGCTCATCAATGTTTACCCAAACACTTCTCAATATGGACCAAAAGATC
ACCAAGTGTAAACCTCTCTGCAATCGAGGCTGGACCACCTGAGTCAGCAGAATCTACAAATATACTCA
AACTTTGTCTCGTGAAGAATTCCTGAGACTGTGTAAAAAAATCATGATGAGATCTATCCAATAAAAAA
GAGAGAGGACCGCAGACGCCTGGCTCTCATCATATGCAATACAAAGTTTGATCACCTGCCTGCAAGGAAT
GGGGCTCACTATGACATCGTGGGGATGAAAAGGCTGCTTCAAGGCCTGGGCTACACTGTGGTTGACGAAA
AGAATCTCACAGCCAGGGATATGGAGTCAGTGTGAGGGCATTGCTGCCAGACCAGAGCACAAGTCCTC
TGACAGCACGTTCTTGGTACTCATGTCTCATGGCATCTAGAGGGAATCTGCGGAACCTGCGCATAAAAAAG
AAAAAACCGGATGTGCTGCTTTATGACACCATCTTCCAGATATTCAACAACCGCAACTGCCTCAGTCTAA
AGGACAAACCAAGGTATCATTGTCCAGGCCTGCAGAGGTGAAAAACATGGGGAACCTCTGGGTCAGAGA
CTCTCCAGCATCTTGGCACTCATCTCTTCCAGTCATCTGAGAACCTGGAGGCAGATTCTGTTTGAAG
ATCCACGAGGAGAAGGACTTCATTGCTTTCTGTTCTTCAACACCACATAACGTGTCTGGAGAGACCGCA
CAAGGGGCTCCATCTTACGGAACCTCATCATGCTTCCAGAAATATTCTTGCTGCTGCCACCTAAT
GGAAATATTTGGAAGGTACAGAAATCATTTGAAGTTCCACAGGCTAAAGCCAGATGCCACCATAGAA
CGAGCAACCTTGACAAGAGATTTCTACCTCTTCTGGAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>RC229221 representing NM_004347 Red=Cloning site Green=Tags(s) MAEDSGKKRRKNFEAMFKGILQSGLDNFVINHMLKNNVAGQTSIQTLVPNTDQKSTSVKKDNHKKKTVK MLEYLKGDV LHGVFN YLAKHDVLT LKEEEKKKYYDTKIEDKAL ILVDSL RKNRVAHQMFQTLLNMDQKI TSVKPLLQIEAGPPESAESTN ILKLCPREEF LRLCKKNHDEIYPIKKREDRRRLALIICNTKFDHLPARN GAHYDI VGMKRLLQGLGYTVVDEKNLTARDMESV LRAFAARPEHKSSDSTFLVLM SHGILEGICGTAHKK KKPDVLLYDTIFQIFNNRNCLSLKDKPKV IIVQACRGEKHGELWVRDSPASLALISSQSENLEADSVCK IHEEKDFIAFCSTPHNVSWDRTRGSIFITELITCFQKYSCCCHLMEIFRKVQKSFEVPQAKAQMP TIE RATLTRDFYLFPGN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004347
ORF Size:	1302 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004347.3
RefSeq ORF:	1305 bp
Locus ID:	838
UniProt ID:	P51878
Cytogenetics:	11q22.3
MW:	49.6 kDa