

Product datasheet for **RC229276**

ALDH7A1 (NM_001182) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALDH7A1 (NM_001182) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ALDH7A1
Synonyms:	ATQ1; EPD; PDE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC229276 representing NM_001182
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTGGCGCCTTCTCGCGCGTGTGTGTGCACGCTGCAAGACCAGCAAGCTCTCTGGACCTTGGAGCA
GGCCTGCCGCTTCATGTCCACTCTCTCATCAATCAGCCCCAGTATGCGTGGCTGAAAGAGCTGGGGCT
CCGCGAGGAAAACGAGGGCGTATAATGGAAGCTGGGGAGGCCGGGAGAGGTTATTACGACCTATTGC
CCTGCTAAACAACGAGCCAATAGCAAGAGTCCGACAGGCCAGTGTGGCAGACTATGAAGAACTGTAAAGA
AAGCAAGAGAAGCATGGAATCTGGGCAGATATCTCTGCTCCAAAACGAGGAGAAATAGTAAGACAGAT
TGGCGATGCCTTGGGGAGAAGTCCAAGTACTAGGAAGCTTGGTGTCTTTGGAGATGGGAAAACTTA
GTGGAAGTGTGGGTGAAGTTCAGGAGTATGTGGATATCTGTGACTATGCTGTTGGTTTATCAAGGATGA
TTGGAGGACCTATCTTGCCTTCTGAAAGATCTGGCCATGCACTGATTGAGCAGTGAATCCCGTAGGCCCT
GGTTGGAATCATCAGGCATTCAATTTCCCTGTGGCAGTGTATGGTTGGAACAACGCCATGCCATGATC
TGTGGAATGTCTGCCTCTGAAAGGAGCTCCAACCCTCCCTCATTAGTGTGGCTGTCAAAAAGATAA
TAGCCAAGTTCTGGAGGACAACAAGCTGCCTGGTGAATTTGTTCTTACTTGTGGTGGAGCAGATAT
TGGCACAGCAATGGCCAAAGATGAACGAGTGAACCTGCTGTCTTCACTGGGAGCACTCAGGTGGAAAA
CAGGTGGGCTGATGGTGCAGGAGAGGTTTGGGAGAAGTCTGTTGGAAGTGGAGGAAACAATGCCATTA
TTGCCTTTGAAGATGCAGACCTCAGCTTAGTTGTTCCATCAGCTCTTTCGCTGCTGTGGGAACAGCTGG
CCAGAGGTGTACCACTGCGAGGCGACTGTTTATACATGAAAGCATCCATGATGAGGTTGTAACAGACTT
AAAAGGCCTATGCACAGATCCGAGTTGGGAACCCATGGGACCCTAATGTTCTCTATGGCCACTCCACA
CCAAGCAGCAGTGAGCATGTTTCTTGGAGCAGTGGGAAGAACAAAGAAAGAGTGGCAGTGGTCTA
TGGGGCAAGGTTATGGATCGCCCTGAAATTATGTAGAACCACAATTGTGACAGGTCTTGGCCACGAT
CGTCCATTGCACACACAGAGACTTTTGTCTCCGATTCTCTATGTCTTTAAATTCAAGAATGAAGAAGAGG
TCTTTGCATGGAATAATGAAGTAAAACAGGACTTTCAAGTAGCATCTTTACAAAGATCTGGCAGAAT
CTTTTCGCTGGCTTGGACCTAAAGGATCAGACTGTGGCATTGTAATGTCAACATTCACAAAGTGGGGCT
GAGATTGGAGGTGCCTTTGGAGGAGAAAAGCACACTGGTGGTGGCAGGGAGTCTGGCAGTGTGCTGGA
AACAGTACATGAGAAGGTCTACTTGTACTATCAACTACAGTAAAGACCTTCTCTGGCCCAAGGAATCAA
GTTTCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC229276 representing NM_001182
Red=Cloning site Green=Tags(s)

MWRLPRALCVHAAKTSKLSGPWSRPAAFMSTLLINQPQYAWLKELGLRENEGVYNGSWGGRGEVITTYC
PANNEPIARVRQASVADYEETVKKAREAWKIWADIPAPKRGEIVRQIGDALREKIQVLGSLVLEMGKIL
VEGVGEVQEYVDICDYAVGLSRMIGGPIILPERSGHALIEQWNPVGLVGIITAFNFPVAVYGWNNAIAMI
CGNVCLWKGAPTTSLISVAVTKIIIAKVLEDNKLPGAICSLTCGGADIGTAMAKDERVNLSTFTGSTQVQK
QVGLMVQERFGRSLELGGNNAIIAFEDADLSLVVPSALFAAVGTAGQRCTTARRLFIHESIHDVNNRL
KKAYAQIRVGNPWPVNLVYGPLHTKQAVSMFLGAVEEAKKEGGTVVYGGKVMRPGNYVEPTIVTGLGHD
ASIAHTETFAPILYVFKFKNEEEVFAWNNEVKQGLSSSIFTKDLGRIFRWLGPKGSDCGIVNVIPTSGA
EIGGAFGGKHTGGRESGSDAWKQYMRRSTCTINYSKDLPLAQGIKFK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8066_c04.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:

ACCN: NM_001182

ORF Size: 1617 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001182.5](#)

RefSeq ORF:	1620 bp
Locus ID:	501
UniProt ID:	P49419
Cytogenetics:	5q23.2
Domains:	aldedh
Protein Families:	Druggable Genome
Protein Pathways:	Arginine and proline metabolism, Ascorbate and aldarate metabolism, beta-Alanine metabolism, Butanoate metabolism, Fatty acid metabolism, Glycerolipid metabolism, Glycolysis / Gluconeogenesis, Histidine metabolism, Limonene and pinene degradation, Lysine degradation, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism, Tryptophan metabolism, Valine, leucine and isoleucine degradation
MW:	58.3 kDa
Gene Summary:	The protein encoded by this gene is a member of subfamily 7 in the aldehyde dehydrogenase gene family. These enzymes are thought to play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. This particular member has homology to a previously described protein from the green garden pea, the 26g pea turgor protein. It is also involved in lysine catabolism that is known to occur in the mitochondrial matrix. Recent reports show that this protein is found both in the cytosol and the mitochondria, and the two forms likely arise from the use of alternative translation initiation sites. An additional variant encoding a different isoform has also been found for this gene. Mutations in this gene are associated with pyridoxine-dependent epilepsy. Several related pseudogenes have also been identified. [provided by RefSeq, Jan 2011]

Product images:



Circular map for RC229276



Western blot validation of overexpression lysate (Cat# [LY432295]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC229276 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).