

## Product datasheet for **RC229284L3V**

### ZNF415 (NM\_001136038) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF415 (NM_001136038) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF415
Synonyms:	Pact; ZfLp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001136038
ORF Size:	1665 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229284).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001136038.2</a> , <a href="#">NP_001129510.2</a>
RefSeq Size:	2547 bp
RefSeq ORF:	1668 bp
Locus ID:	55786
UniProt ID:	<a href="#">Q09FC8</a>
Cytogenetics:	19q13.42
Protein Families:	Transcription Factors
MW:	64 kDa



[View online »](#)

**Gene Summary:**

Involved in transcriptional regulation. Transcriptional activity differed among the various isoforms. All isoforms except isoform 3 seem to suppresses the transcriptional activities of AP-1 and p53/TP53.[UniProtKB/Swiss-Prot Function]