

Product datasheet for **RC229315**

DACT3 (NM_145056) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DACT3 (NM_145056) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DACT3
Synonyms:	DAPPER3; RRR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC229315 representing NM_145056.
 Blue=ORF Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCCGGGCCCTTCTCGTTCCCCGGTGAGCCCTGAGCGGGCCGGCTGCGGGGCTGGCTGGAGGGTAGC
CTGGCCGGGCTCTGCGAGTTACATTGGCTCCGGGAGAGCAGGAGTACCGCGTGCAGCAGCGCTGCGG
CTGGCCAGCCCGAATGGGGGGCGCCGAGGCCGAGGACGAGGAGACCGCGATGAGGATGAAGATGCG
GCGGCGCGCGCCGGCCGACGCGCCCTGGAGGAGCAGCTGGAGGCCCTGCCTGGTCTCGTCTGGGAC
CTGGGACAGCAGCTGGGAGACCTGAGCCTGGAGTCTGGGGCCTGGAACAGGAGAGCGGGCGTAGCTCG
GGCTTCTATGAAGATCCAGCTCTACAGGAGGTCCAGATTACCACCCTCAACCTTCTGTGGGACAGT
GGCTTCTCTGGATCCAGCTCCTATGGTCGCTGGGTCCCTCTGAGCCCCGGGCATCTATGCCAGTGAG
AGGCCAAGTCCCTAGGAGACGCCAGTCCAGCGCTCCGAGGTGGTGGCGCGCGGGCAGCGGTGCCG
CGGTCTTCTCAGCGCCCTACCCGACGGCAGGGGGTCCGCCGGCCGGAGGCCCTGCTCCTCGCGGAG
CGCGGGGCCCGCCGGGCCCTTCTGACGCCAGCCCCCTGCAGCCGTGGCGATGCGCAGCCCGCGG
CCCTGCGGCCGCCCTCCACCGACTCGCCGACGCGGGGGCGCAGGGCGGCCCTGGACGGCTACATC
TCGGCGCTCCTGCGCAGGCGCCCGCCGGGGGGCGGGCCAGCCCCGGACCAGTCCCGGGGGCGCGGAC
GGCGGCCCGCGCGCCAGAACAGCGTGCGCCAGCGGCCGCCGACCGCTCTCCGTCCCCCGGACGCGG
CGACCCGCGCGGGAGCCCTCGTTGGAGCGCTCGGGGGCCACCCACCAGCCCTGCCGCCTTGAGCCGC
GCCTGGCGTCTGCTGGGAGTCGGAGGCGGCACCCGAGCCCGCTGCGCCGCCCGCCGCCCTCACCC
CCCAGACCCCGCTGAGGGCCGCTTGGTGAAGGCGCAGTACATCCGGGCGCGCAGGCGGCCACCCGA
GGCTCCTTGGCCGCGCCCGCCGCCAAACCCCGCCACTGACCCGCGCCGACGCTGGAGCAGTCA
CACCCCGGAGCGTCCCGGGCCCGCGCCCGCTGGACGCATGGCCGAGGCTTCGGGGCCCGCGGC
TCGCCCAGGGCCCGCAAGGCCCTCGCGCTCCAGTCTGAGACCAGCCTGCTGGGCCCGCCCTCCGCGGTC
CCTTCGGGGCCCTAAGTACCCACGCGGAGCGGGAAGAGCCTCGGCCCTCACGGCCACGCGCGGC
CCAGCGCCACGCTGGCGGCCAGGCCGCGAGGTCCTGCCGTGCTGGCGCTCCACTGCGGAGATCGAC
GCTGCCGATGGGCGCGCGTGGCCCGGAGCCCTGCGGCGCGTGTCCCGGCCCGCGCCGTCGCCG
TCAGCTCCCAGCGTCTGCTTTACGGCTGCGCGGCAGCGACTCCGAGTCTCGGCTGGGCGCTG
GGGCCCTGGGACCGCGGGGCCCTGCGGGAGCGCTCGGCGGGGTTACGGGAGAGCGAATCGAGCGCC
AGCGAGGAGAATCGCTGCCTCAGCTCTGCCTCCAGCGACTCAGACGGCAGCGGTGCCCTCGTGTGG
CCGACGAGCTGGTGGCGCCACCGCGCCTCTGGGGTGGAGCAGGTGCAGGGCGCCCGCAGGCCCC
GCCAAAGTCTTCGTGAAATCAAAGCTCCACGCGCTCAAGAAAAGATACTGCGTTTCGGTTCGGGT
TCTCTAAGGTCATGACTACAGTGACGGACCGACGCGTACGCGGCCCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTAA
  
```

Protein Sequence:

>Peptide sequence encoded by RC229315
 Blue=ORF Red=Cloning site Green=Tag(s)

```


MIRAFSFPVSPERGLRGWLEGLAGLCELHWRERQYRVQALRLAQPGMGGAEEDEEDA
AAARRAAALEEQLEALPGLVWDLGQQLGDL SLESGLEQESGRSSGFYEDPSSTGGPDSPPSTFCGDS
GFSGSSSYGRLGPSEPRGIYASERPKSLGDASPSAPEVVGARAAVPRSFSAFYPTAGGSAGPEACSSAE
RRARAGPFLTPSPLHAVAMRSRPRCGRPPTDSPDAGGAGRPLDGYISALLRRRRRRGAGQPRTPGGAD
GGPRRQNSVRQPPDASPSGARSAPREPSLERVGGHPTSPAALSRAWASSWESEAPEAAPPAAAPSP
PDSPAEGRLVKAQYIPGAQAATRGLPGRAARRKPPPLTRGRSVEQSPPRERPRAGRRRMAEASGRRG
SPRARKASRSQSETSLLGRASAVPSGPPKYPTAEREPRPRRGPAPTLAAQAAGSCRWRSTAEID
AADGRRVRPRAPAARVPGPGSPSAPQRLLYGCAGSDSECSAGRLGPLGRRGPAGGVGGYGESESSA
SEGESPAFSSASSSDSGGLVWPQQLVAATAASGGGAGAGAPAGPAKVFVFIKASHLKKKILFRFRSG
SLKVMTTV
SGPTRRTRPLEQKLISEEDLAANDILDYKDDDDKV
  
```

Chromatograms:

https://cdn.origene.com/chromatograms/mk8070_e09.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_145056

ORF Size: 1887 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145056.3](#)

RefSeq ORF: 1890 bp

Locus ID: 147906


UniProt ID: [Q96B18](#)

Cytogenetics: 19q13.32

MW: 64.9 kDa

Gene Summary: May be involved in regulation of intracellular signaling pathways during development. Specifically thought to play a role in canonical and/or non-canonical Wnt signaling pathways through interaction with DSH (Dishevelled) family proteins.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC229315