

## Product datasheet for RC229322L4V

### SLC2A13 (NM\_052885) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC2A13 (NM_052885) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC2A13
Synonyms:	HMIT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052885
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229322).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_052885.3</a>
RefSeq ORF:	1947 bp
Locus ID:	114134
UniProt ID:	<a href="#">Q96QE2</a>
Cytogenetics:	12q12
Domains:	sugar_tr
Protein Families:	Transmembrane
MW:	70.2 kDa


[View online »](#)

**Gene Summary:** H(+)-myo-inositol cotransporter (PubMed:11500374). Can also transport related stereoisomers (PubMed:11500374).[UniProtKB/Swiss-Prot Function]