

Product datasheet for **RC229364L3V**

ZYG11A (NM_001004339) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZYG11A
Synonyms:	ZYG11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001004339
ORF Size:	2277 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC229364).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001004339.2
RefSeq ORF:	2280 bp
Locus ID:	440590
UniProt ID:	Q6WRX3
Cytogenetics:	1p32.3
MW:	85.6 kDa



Gene Summary:

Probably acts as target recruitment subunit in an E3 ubiquitin ligase complex ZYGA-CUL2-elongin BC.[UniProtKB/Swiss-Prot Function]