

## Product datasheet for **RC229383L3V**

### ACAP3 (NM\_030649) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ACAP3 (NM_030649) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACAP3
Synonyms:	CENTB5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030649
ORF Size:	2502 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229383).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030649.2</a>
RefSeq ORF:	2505 bp
Locus ID:	116983
UniProt ID:	<a href="#">Q96P50</a>
Cytogenetics:	1p36.33
Domains:	ArfGap, PH, ANK
Protein Pathways:	Endocytosis
MW:	92.3 kDa



[View online »](#)

**Gene Summary:** GTPase-activating protein for the ADP ribosylation factor family.[UniProtKB/Swiss-Prot Function]