

Product datasheet for **RC229415L3V**

ZCCHC2 (NM_017742) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZCCHC2
Synonyms:	C18orf49
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017742
ORF Size:	3534 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC229415).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_017742.4 , NP_060212.4
RefSeq Size:	5559 bp
RefSeq ORF:	3537 bp
Locus ID:	54877
UniProt ID:	Q9C0B9
Cytogenetics:	18q21.33



Domains: zf-CCHC

MW: 125.9 kDa