

## Product datasheet for **RC229438L2V**

### Diffuse panbronchiolitis critical region protein 1 (DPCR1) (NM\_080870) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Diffuse panbronchiolitis critical region protein 1 (DPCR1) (NM_080870) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MUCL3
Synonyms:	C6orf37; DPCR1; PBLT
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_080870
ORF Size:	4179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229438).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080870.3</a> , <a href="#">NP_543146.2</a>
RefSeq ORF:	4182 bp
Locus ID:	135656
UniProt ID:	<a href="#">Q3MIW9</a>
Cytogenetics:	6p21.33
Protein Families:	Transmembrane
MW:	151 kDa



[View online »](#)

**Gene Summary:**

May modulate NF-kappaB signaling and play a role in cell growth.[UniProtKB/Swiss-Prot  
Function]