

## Product datasheet for **RC229439L3V**

### SCAPER (NM\_020843) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SCAPER (NM_020843) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCAPER
Synonyms:	IDDRP; MSTP063; Zfp291; ZNF291
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020843
ORF Size:	4200 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229439).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020843.2</a>
RefSeq ORF:	4203 bp
Locus ID:	49855
UniProt ID:	<a href="#">Q9BY12</a>
Cytogenetics:	15q24.3
Domains:	ZnF_U1
MW:	158.1 kDa


[View online »](#)

**Gene Summary:** CCNA2/CDK2 regulatory protein that transiently maintains CCNA2 in the cytoplasm.  
[UniProtKB/Swiss-Prot Function]