

## Product datasheet for **RC229499A1V**

### Human **NDUFC1 (NM\_001184986)** AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human NDUFC1 (NM\_001184986) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** NDUFC1  
**Synonyms:** KFYI  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC229499 representing NM\_001184986  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

**ATGGCGCCGTCCGCCTTGCTGCGTCCCCTTTCCCGGCTGCTGGCCCCGCCAGGCTCCCGAGCGGCCCTT**  
**CAGTGCGATCAAAGTTCTACGTGCGAGAGCCCGCAATGCCAAACCTGACTGGCTGAAAGTTGGGTTAC**  
**CTTGGGCACTGTCTTCTTGGATCTATCTCATCAAACAACAATGAAGATATTTAGAGTACAAA**  
**AGAAGAAATGGGCTGGAA**

**ACGCGTACGCGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT**  
**ACAAGGATGACGACGATAAGGTTTAA**

**Protein Sequence:** >RC229499 representing NM\_001184986  
**Red**=Cloning site **Green**=Tags(s)

MAPSALLRPLSRLAPARLPSGPSVRSKFYVREPPNAKPDWLKVGFTLGTTVFLWIYLIKQHNEDILEYK  
RRNGLE

**TRTRPLEQKLI SEEDLAANDILDYKDDDDKV**

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001184986  
**ORF Size:** 228 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

**Stability:** AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

**RefSeq:** [NM\\_001184986.1](#), [NP\\_001171915.1](#)

**RefSeq Size:** 1257 bp

**RefSeq ORF:** 231 bp

**Locus ID:** 4717

**UniProt ID:** [O43677](#)

**Cytogenetics:** 4q31.1

**MW:** 8.7 kDa