

## Product datasheet for **RC229505L3V**

### NDUFC1 (NM\_001184991) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NDUFC1 (NM_001184991) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDUFC1
Synonyms:	KFYI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001184991
ORF Size:	228 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229505).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001184991.1</a> , <a href="#">NP_001171920.1</a>
RefSeq Size:	1198 bp
RefSeq ORF:	231 bp
Locus ID:	4717
UniProt ID:	<a href="#">O43677</a>
Cytogenetics:	4q31.1
Protein Pathways:	Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease



[View online »](#)

**MW:** 8.7 kDa

**Gene Summary:** The encoded protein is a subunit of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010]