

## Product datasheet for RC229512

### TRIQK (NM\_001171795) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TRIQK (NM\_001171795) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TRIQK  
**Synonyms:** C8orf83; PRO0845; UPF0599  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC229512 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGTAGAAAAGATGCTGCTACTATAAACTTCTGTTGATCAGTACAGAAAACAAATTGGTAAACAGG  
ATTATAAAAAAATAAACCTATTTACGAGCAACCAAATTAAGCAGAAGCAAAGAAAACAGCAATAGG  
CATAAAGGAAGTTGGCCTTGTACTTGCAGCTATATTGGCACTACTACTGGCTTTCTATGCTTTCTTTAT  
CTCAGACTCACCACGGATGTTGACCCTGATCTGGACCAAGATGAAGAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC229512 protein sequence  
Red=Cloning site Green=Tags(s)  
MGRKDAATIKLPVDQYRKQIGKQDYKTKPILRATKCLKAEAKKTAIGIKEVGLVLAAILALLLAFYAFFY  
LRLTTDVPDLQDED

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8036\\_h12.zip](https://cdn.origene.com/chromatograms/mk8036_h12.zip)

**Restriction Sites:** Sgfl-MluI



[View online »](#)

## Cloning Scheme:



ACCN: NM\_001171795

ORF Size: 258 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001171795.2](#)


RefSeq Size: 3684 bp

RefSeq ORF: 261 bp

Locus ID: 286144

UniProt ID: [Q629K1](#)  
 Cytogenetics: 8q22.1  
 MW: 9.7 kDa  
 Gene Summary: May play a role in cell growth and maintenance of cell morphology.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RC229512