

Product datasheet for **RC229512A1V**

Human TRIQK (NM_001171795) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TRIQK (NM_001171795) AAV Particle
Tag: Myc-DDK
Symbol: TRIQK
Synonyms: C8orf83; PRO0845; UPF0599
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC229512 representing NM_001171795
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGTAGAAAAGATGCTGCTACTATAAACTCCTGTTGATCAGTACAGAAAACAAATTGGTAAACAGG
ATTATAAAAAAATAACCTATTTTACGAGCAACCAATTTAAAGCAGAAGCAAAGAAAACAGCAATAGG
CATAAAGGAAGTTGGCCTTGTACTTGCAGCTATATTGGCACTACTGCTTTCTATGCTTTCTTTTAT
CTCAGACTCACCAAGGATGTTGACCCTGATCTGGACCAAGATGAAGAT

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229512 representing NM_001171795
Red=Cloning site **Green**=Tags(s)

MGRKDAATIKLPVDQYRKQIGKQDYKTKPILRATKLAEAKKTAIGIKEVGLVLAAILALLLAFYAFFY
LRLTTDVPDLDDQDED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001171795
ORF Size: 258 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001171795.1</u>
RefSeq Size:	3684 bp
RefSeq ORF:	261 bp
Locus ID:	286144
UniProt ID:	<u>Q629K1</u>
Cytogenetics:	8q22.1
MW:	9.7 kDa