

Product datasheet for **RC229515A1V**

Human BLCAP (NM_001167821) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human BLCAP (NM_001167821) AAV Particle
Tag: Myc-DDK
Symbol: BLCAP
Synonyms: BC10
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC229515 representing NM_001167821
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTATTGCCTCCAGTGGCTGCTGCCCGTCTCCTCATCCCCAAGCCCCTCAACCCCGCCCTGTGGTTCA
GCCACTCCATGTTTCATGGGCTTCTACCTGCTCAGCTTCTCCTGGAACGGAAGCCTTGACAATTTGTGC
CTTGGTTTTCTGGCAGCCCTGTTCTTATCTGCTATAGCTGCTGGGGAAACTGTTTCTGTACCACTGC
TCCGATCCCCGCTTCCAGAATCGGCGCATGATCCCGCGTTGTGGGCACC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229515 representing NM_001167821
Red=Cloning site Green=Tags(s)

MYCLQWLLPVLLIPKPLNPALWFSHSMFMGFYLLSFLLEKPKPTICALVFLAALFLICYSCWGNCFLYHC
SDSPLPESAHDPGVVGT

TRTRPLEQKLI**SEEDLA**ANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001167821
ORF Size: 261 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001167821.1</u> , <u>NP_001161293.1</u>
RefSeq Size:	2046 bp
RefSeq ORF:	264 bp
Locus ID:	10904
UniProt ID:	<u>P62952</u>
Cytogenetics:	20q11.23
MW:	9.9 kDa