

## Product datasheet for **RC229517L4V**

### **BLCAP (NM\_001167823) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | BLCAP (NM_001167823) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | BLCAP  |
| Synonyms:                 | BC10   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001167823   |
| ORF Size:                 | 261 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC229517).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001167823.1</a>   |
| RefSeq Size:              | 2119 bp  |
| RefSeq ORF:               | 264 bp   |
| Locus ID:                 | 10904  |
| UniProt ID:               | <a href="#">P62952</a>   |
| Cytogenetics:             | 20q11.23   |
| Protein Families:         | Transmembrane  |
| MW:                       | 9.9 kDa  |



[View online »](#)

**Gene Summary:**

This gene encodes a protein that reduces cell growth by stimulating apoptosis. Alternative splicing and the use of alternative promoters result in multiple transcript variants encoding the same protein. This gene is imprinted in brain where different transcript variants are expressed from each parental allele. Transcript variants initiating from the upstream promoter are expressed preferentially from the maternal allele, while transcript variants initiating downstream of the interspersed NNAT gene (GeneID:4826) are expressed from the paternal allele. Transcripts at this locus may also undergo A to I editing, resulting in amino acid changes at three positions in the N-terminus of the protein. [provided by RefSeq, Nov 2015]