

## Product datasheet for **RC229554A1V**

### Human Serum Amyloid A (SAA1) (NM\_001178006) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human Serum Amyloid A (SAA1) (NM\_001178006) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Serum Amyloid A  
**Synonyms:** PIG4; SAA; SAA2; TP53I4  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC229554 representing NM\_001178006  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAAGCTTCTCACGGGCTGGTTTTCTGCTCCTTGGTCTGGGTGTCAGCAGCCGAAGCTTCTTTTCGT  
TCCTTGGCGAGGCTTTTGATGGGGCTCGGGACATGTGGAGAGCCTACTCTGACATGAGAGAAGCCAATTA  
CATCGGCTCAGACAAATACTTCCATGCTCGGGGAACATGATGCTGCCAAAAGGGGACCTGGGGGTGTC  
TGGGCTGCAGAAGCGATCAGCGATGCCAGAGAGAATATCCAGAGATTCTTTGGCCATGGTGCGGAGGACT  
CGCTGGCTGATCAGGTCCTCAATGAATGGGGCAGGAGTGGCAAAGACCCCAATCACTTCCGACCTGCTGG  
CCTGCCTGAGAAATAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC229554 representing NM\_001178006  
Red=Cloning site Green=Tags(s)

MKLLTGLVFCSLVLGVSSRSFFSFLGEAFDGDARMWRAYSMDREANYIGSDKYFHARGNYDAAKRGPGGV  
WAAEAISDARENIQRFFGHGAEDSLADQAANEWGRSGKDPNHFRPAGLPEKY

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001178006



[View online >](#)

<b>ORF Size:</b>	366 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001178006.1</a></u>
<b>RefSeq Size:</b>	592 bp
<b>RefSeq ORF:</b>	369 bp
<b>Locus ID:</b>	6288
<b>UniProt ID:</b>	<u><a href="#">P02735</a></u>
<b>Cytogenetics:</b>	11p15.1
<b>MW:</b>	13.5 kDa