

Product datasheet for **RC229555**

FAM229A (NM_001167676) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: Myc-DDK

Symbol: FAM229A

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC229555 representing NM_001167676
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCTGCCCTCTCGACGCCCGGGCCGGGCACGCCACAGAGACCTGCCCGGCTCCGCCTGGACCGGAGC
GTTCTCCCGCGCCAGGGCTCCGGCAGTGTCTTAGCCTGGACCGGTCTCAACCGCCGGGAGAGCGCC
CCGGGGCTTGGACATGAGCGCCAGGAGCCCCGCAGGGTCGGAGATCCCCATTGAGGCCGGAGACTCC
CGTGGCCTTGGCGCGCCCCGAGTCCCAGGACAGCCCGGAGGCGGTAGCGACGGAGACAACCCGGTCA
GGCCGCTTCGTGCTGCCCTGGATGCCACTGCCTGACGCTGCTGCACGTGCCCATCGACGTCTACCTCGC
CATGGCGGGAGCCCCCGGGCCCGCCACG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229555 representing NM_001167676
Red=Cloning site Green=Tags(s)

MLPSSTPGPGHATETCPAPGPERSPAARAPAAASSLGPVSTAGRAPRLDMSAQEPQGRRFPIEAGDS
RGLAAPESQDSPEAVATEHNPVRLRRCPGCHCLTLLHVPIDVYLAMGGSPRARAT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8036_g08.zip

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001167676

ORF Size: 381 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001167676.1, NP_001161148.1
RefSeq ORF: 384 bp
Locus ID: 100128071
UniProt ID: H3BQW9
Cytogenetics: 1p35.1
MW: 13.4 kDa

Product images:



Circular map for RC229555



Coomassie blue staining of purified FAM229A protein (Cat# [TP329555]). The protein was produced from HEK293T cells transfected with FAM229A cDNA clone (Cat# RC229555) using MegaTran 2.0 (Cat# [TT210002]).