

Product datasheet for RC229557

MCFD2 (NM_001171511) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MCFD2 (NM_001171511) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: MCFD2
Synonyms: F5F8D; F5F8D2; LMAN1IP; SDNSF
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC229557 representing NM_001171511
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTCAGCGTCTGCAGCTGCCGGACCAGCTCCGGCATGCGGTCCAGTGGCCCTCGGCGCGGCAGCGCT
CCAGCTCGCTCTCCACCTTCAGGCATATCATGGAGCATCTAGAAGGTGCATCAACAAACCAGAGGCGGA
GATGTCGCCACAAGAATTGCAGTCCATTACTTCAAATGCATGATTATGATGGCAATAATTTGCTTGAT
GGCTTAGAACTCTCCACAGCCATCACTCATGTCCATAAGGAGGAAGGGAGTGAACAGGCACCACTAATGA
GTGAAGATGAACTGATTAACATAATAGATGGTGTTTTGAGAGATGATGACAAGAACAATGATGGATACAT
TGACTATGCTGAATTTGCAAAATCACTGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229557 representing NM_001171511
Red=Cloning site Green=Tags(s)

MLSVCSRTSSGMRSQWPSARQRSSSLSTFRHIMEHLEGVINKPEAEMSPQELQLHYFKMHDYDGNLLD
GLELSTAIHVHKEEGSEQAPLMSDELINIIDGVLRDDDKNNDGYIDYAEFAKSLQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001171511

ORF Size: 381 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001171511.2 , NP_001164982.1
RefSeq Size:	4179 bp
RefSeq ORF:	384 bp
Locus ID:	90411
UniProt ID:	Q8NI22
Cytogenetics:	2p21
MW:	14.4 kDa
Gene Summary:	<p>This gene encodes a soluble luminal protein with two calmodulin-like EF-hand motifs at its C-terminus. This protein forms a complex with LMAN1 (lectin mannose binding protein 1; also known as ERGIC-53) that facilitates the transport of coagulation factors V (FV) and VIII (FVIII) from the endoplasmic reticulum to the Golgi apparatus via an endoplasmic reticulum Golgi intermediate compartment (ERGIC). Mutations in this gene cause combined deficiency of FV and FVIII (F5F8D); a rare autosomal recessive bleeding disorder characterized by mild to moderate bleeding and coordinate reduction in plasma FV and FVIII levels. This protein has also been shown to maintain stem cell potential in adult central nervous system and is a marker for testicular germ cell tumors. The 3' UTR of this gene contains a transposon-like human repeat element named 'THE 1'. A processed RNA pseudogene of this gene is on chromosome 6p22.1. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Apr 2016]</p>