

## Product datasheet for **RC229563A1V**

### Human IFITM10 (NM\_001170820) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human IFITM10 (NM\_001170820) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** IFITM10  
**Synonyms:** DSPA3  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC229563 representing NM\_001170820  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGCCACACTGTTCCCATGGAGTCCAAGAGCAGCAAGACCGACAGCGTGC GGCGTCCGGCGCGC  
CCCCTGCCTGCAAGCACCTAGCCGAGAAGAAGACGATGACCAACCCACGACCGTCATCGAGGTCTACCC  
GGACACCACCGAGGTGAACGACTATTACCTGTGGTCCATCTTCAACTTCGCTACCTCAACTTCTGCTGC  
CTGGGCTTCATCGCCTTGGCCTACTCCCTCAAAGTCCGAGACAAGAAGCTTCTCAATGACCTGAATGGAG  
CCGTGGAGGATGCAAAGACGGCCCGGTGTTCAACATCACCAGTTCTGCCCTGGCAGCCTCCTGCATCAT  
CCTCGTCTTCATCTTCTGCGGTACCCCTCACCGACTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC229563 representing NM\_001170820  
Red=Cloning site Green=Tags(s)

MAPTLFPMESKSSKTDSVRAAGAPPACKHLAEKKTMTNPTTVIEVYPTTTEVNDYYLWSIFNFVYLNFCF  
LGFIALAYSLKVRDKLLNDLNGAVEDAKTARLFNITSSALAASCIILVFIFLRYPLTDY

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001170820



[View online »](#)

<b>ORF Size:</b>	390 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001170820.2</a></u> , <u><a href="#">NP_001164291.1</a></u>
<b>RefSeq ORF:</b>	687 bp
<b>Locus ID:</b>	402778
<b>UniProt ID:</b>	<u><a href="#">A6NMD0</a></u>
<b>Cytogenetics:</b>	11p15.5
<b>MW:</b>	14.52 kDa