

Product datasheet for RC229599L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Aquaporin 1 (AQP1) (NM_001185062) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Aquaporin 1 (AQP1) (NM_001185062) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Aquaporin 1

Synonyms: AQP-CHIP; CHIP28; CO

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC229599).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:






 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.





Plasmid Map:



ACCN: NM_001185062

ORF Size: 462 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001185062.1</u>, <u>NP 001171991.1</u>

RefSeq ORF: 465 bp

Locus ID: 358



Aquaporin 1 (AQP1) (NM_001185062) Human Tagged Lenti ORF Clone - RC229599L4

Cytogenetics: 7p14.3

Protein Families: Druggable Genome, Ion Channels: Other, Transmembrane

MW: 17.1 kDa

Gene Summary: This gene encodes a small integral membrane protein with six bilayer spanning domains that

functions as a water channel protein. This protein permits passive transport of water along an osmotic gradient. This gene is a possible candidate for disorders involving imbalance in

ocular fluid movement. [provided by RefSeq, Aug 2016]