

Product datasheet for **RC229606**

EOLA1 (NM_001171908) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: EOLA1 (NM_001171908) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: EOLA1
Synonyms: CXorf40; CXorf40A
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC229606 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGTTTGGCTGCCTCTCCTCCGGCAGCCTTATGCTGGCTTTGTCTTAAATGGAATCAAGACTGTGG
AGACGCGCTGGCGTCCCCTGCTGAGCAGCCAGCGGAACGTACCATCGCCGTCCACATTGCTCACAGGGA
CTGGGAAGGCGATGCCTGGCGGGAGCTGCTGGTGGAGAGACTCGGGATGACTCCTGCTCAGATTCAAGCC
TTGCTCAGAAAGGGGAAAAGTTTGGTCGAGGAGTGATAGCGGGACTCGTTGACATTGGGAAAACCTTGC
AATGCCCGAAGACTTAACTCCCGATGAGGTTGTGGAAGTAAAGAAATCAAGCTGTACTGACCAACCTGAA
GCAGAAGTACCTGACTGTGATTTCAAACCCAGGTGGTTACTGGAGCCATACCTAGGAAAGGAGGCAAG
GATGTATTCCAGGTAGACATCCCAGAGCACCTGATCCCTTGGGGCATGAAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229606 protein sequence
Red=Cloning site Green=Tags(s)

MKFGCLSFRQPYAGFVNLNGIKTVETRWRLPSSQRNCTIAVHIAHRDWEGDWRELLVERLGMTPAQIQAL
LLRKGEKFRGVIAGLVDIGETLQCPEDLTPDEVLELQAVLTLNKQKYLTVISNPRWLLLEIPRKGK
DVFQVDIPEHLIPLGHEV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6357_a08.zip



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001171908

ORF Size: 474 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171908.3](#)


RefSeq Size: 2404 bp

RefSeq ORF: 477 bp

Locus ID: 91966

UniProt ID: [Q8TE69](#)
 Cytogenetics: Xq28
 MW: 17.9 kDa
 Gene Summary: May have an important role of cell protection in inflammation reaction.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC229606